



All AME Products have CE, ISO 13485, FDA and other quality certificates.

We also have GMP, Free Sales Certificates and complete product registration dossiers.

Solutions to all your
Healthcare Needs

HEALTHCARE CATALOGUE

Our Core Values

At AME, our reach transcends borders, and our commitment to improving global health finds resonance in 17 countries around the world. We appreciate the diversity of our customer base and strive constantly to serve the differentiated needs of each market.


Come explore our expansive range and join us in our journey to redefine the landscape of home-care health solutions across the globe. Trust AME, your Medical Devices Brand, to make health and wellness an attainable goal for all.

Presenting our distinguished manufacturing partners, pivotal in crafting the exceptional quality associated with our AME-branded healthcare products.

“ AME brand is synonymous with premium customer service.

- Dr. M. Tayyab Mustafa, CEO

**AME IS THE LARGEST
HOME CARE
MEDICAL
DEVICES
BRAND
TO ORIGINATE FROM
THE MIDDLE EAST**





Welcome to AME {Ayy-Em-Eee)
Revolutionizing healthcare
through innovative, and cost-
effective healthcare solutions

Our dedication to revolutionizing healthcare has been at the forefront of our operations since our inception in 2004. Born out of a passion for better healthcare, our products thrive at the intersection of quality, efficacy, and affordability; catering to an array of needs with unwavering performance.

Establishing our presence in the Home Care domain, we offer an impressive array of over 200 SKUs spanning across 15 crucial categories. This broad spectrum allows us to address various health and wellness concerns, thus living up to our mission of providing comprehensive healthcare solutions for every home leading to “A better Life”.



EFFICACY

All our products should be effective in delivering the solutions they were manufactured for. We work tirelessly to ensure all our products provide a seamless user experience.



RELIABILITY

We use the highest quality of raw materials procured from the best suppliers. These are processed in state of the art manufacturing sites with staff with decades of experience ensuring that our products are reliable. The product works consistently over time without breakdowns.



AFFORDABILITY

Price always matters. We strive to keep the prices affordable for all our customers. AME also provides segmented price points to cater to all patients.

Our Manufacturing Facility

AME exclusively collaborates with cutting-edge production facilities that serve as the hallmark of excellence. These state-of-the-art factories, where AME Health Care products come to life, embody a dedication to innovative manufacturing. Substantial investments in research and development have yielded a series of patented technologies, placing us at the forefront of advancements in healthcare.

All operations, conducted on an impressive scale, adhere rigorously to quality standards, resulting in the acquisition of multiple quality certificates that underscore our commitment to crafting AME products with unparalleled excellence.

All of our esteemed manufacturing facilities play a pivotal role in the creation of AME-branded healthcare solutions, exemplifying a commitment to quality, innovation, and the provision of cutting-edge products that bear the hallmark of excellence.



Our History




Table of CONTENTS

01



Vital Signs Monitoring Solutions

8

- Blood Pressure Monitors
- Blood Glucose Monitors
- Clinical Essentials
- Digital Thermometers
- Weighing Scales

9
16
18
27
30

02




Respiratory Care Solutions

35

- Nebulizers
- Steam Inhalers
- Humidifiers

36
51
53

03



Rapid Diagnostic Solutions

55

- Pregnancy & Fertility
Test Kits
- Drug of Abuse Kits

56
60

62

Maternity Care Solutions

63

Breast Pumps



04

68

Pain Relief & Wellness Solutions

69
72

Pain & Cooling Patches
Foot Cream



05

73

Safe Injections Solution

74

Auto Disable Syringes



06

75

Surgical Tools & Hospital Furniture Solutions

76
76

Surgical Instruments
Hospital Furniture




07

Vital Signs Monitoring Solutions



At AME, we understand that accurate measurement is the foundation of quality healthcare. Our Vital Signs Monitoring range is designed to deliver precise, reliable results that empower both healthcare professionals and individuals to track essential health indicators with confidence. Our monitoring solutions ensure that critical health data is always accessible; supporting early detection, better decision-making, and improved patient outcomes.

Product Range



Blood Pressure Monitors



Blood Glucose Monitors



Digital Thermometers



Clinical Essentials



Weighing Scales



BLOOD PRESSURE MONITORS

AME Blood Pressure Monitors are designed to empower individuals in managing their health with confidence and ease. Whether you're tracking daily fluctuations or maintaining long-term health records, our AME monitors provide a seamless experience with intuitive interfaces, smart memory functions, and user-friendly designs. With features like large digital displays, multi-user modes, and portable convenience, AME Blood Pressure Monitors ensure health monitoring is both simple and accessible.



6 STAR TECHNOLOGY

Our Digital Blood Pressure Monitors feature advanced 6 Star Smart Technology, ensuring precise measurements, ease of use, and reliable performance every time. Designed for both professionals and home users, this innovation combines accuracy, safety, and comfort for smarter health monitoring.



SMART MOTION

To detect movement during measurement



SMART MEASURE

Reduce measurement time by measuring during inflation



SMART CUFF

To detect if cuff is loose during measurement



SMART DETECT

Irregular heartbeat indicator



SMART INDICATE

3 COLOUR Hypertension Indicator



SMART AFIB

Non Invasive detection of Atrial Fibrillation

BLOOD PRESSURE MONITORS

6 STAR TECHNOLOGY



Monitoring
Solutions



MODEL: AME 1359

FEATURES

- Automatic Power Off
- 2×60 memories with date and time
- Average last 3 results
- Deluxe carry case
- Bluetooth® optional
- Talking optional

DESIGN ELEMENTS

- Silver White Finish
- Blue Backlight
- Large Display
- 4×AAA batteries or AC Adapter Or USB C Type

SPECIFICATIONS

- LCD dimension: 6.7x6.9cm
- Unit dimension: 13.1x10.1x4.4cm

TECHNOLOGY



**SMART
MOTION**



**SMART
MEASURE**



**SMART
CUFF**



**SMART
DETECT**



**SMART
INDICATE**



**SMART
AFIB**



Monitoring
Solutions

BLOOD PRESSURE MONITORS

5 STAR TECHNOLOGY



MODEL: AME 6195



MODEL: AME 6192

TECHNOLOGY



Smart MOTION Smart MEASURE Smart CUFF Smart DETECT Smart INDICATE

FEATURES

- Automatic Power Off
- 2x150 memories with date and time
- Average last 3 results
- Deluxe carry case
- Bluetooth® optional
- Talking optional
- Backlight optional

Low Battery detection

DESIGN ELEMENTS

- Large Display
- 3xAAA batteries or Type-C

4 Colors available:

4 Colors available:

SPECIFICATIONS

LCD dimension: 4.3x6.6cm
Unit dimension: 15.0x8.0x4.6cm

LCD dimension: 5.2x5.2cm
Unit dimension: 14.0x10.3x5.1cm

BLOOD PRESSURE MONITORS

5 STAR TECHNOLOGY



Monitoring
Solutions



MODEL: AME 6196



MODEL: AME 6191

TECHNOLOGY



Smart MOTION Smart MEASURE Smart CUFF Smart DETECT Smart INDICATE

FEATURES

- Automatic Power Off
- 2x150 memories with date and time
- Average last 3 results
- Deluxe carry case
- Bluetooth® optional
- Talking optional
- Backlight optional
- Low Battery detection

DESIGN ELEMENTS

- Large Display
- 3xAAA batteries or Type-C

4 Colors available:

Silver White finish
USB Connector

SPECIFICATIONS

LCD dimension: 6.2x4.7cm
Unit dimension: 10.3x9.3x5.1cm

LCD dimension: 4.6x6.2cm
Unit dimension: 14.0x9.8x4.8cm



Monitoring
Solutions

BLOOD PRESSURE MONITORS

5 STAR TECHNOLOGY



MODEL: AME BP 807

MODEL: AME DBP 6194

MODEL: AME BP 816

TECHNOLOGY



Smart MOTION Smart MEASURE Smart CUFF Smart DETECT Smart INDICATE

FEATURES

Automatic Power Off
Average last 3 results
Deluxe carry case
Bluetooth® optional
Talking optional
Backlight optional

2×60 memories with date and time

2×1500 memories with date and time
Low Battery detection


2×90 memories with date and time
Low Battery detection

DESIGN ELEMENTS

Large Display

4 Colors available: 
4×AAA batteries or AC Adaptor

White finish (colors available with MOQ)
3×AAA batteries or AC Adaptor
USB Connector

4 Colors available: 
4×AAA batteries or Type-C

SPECIFICATIONS

LCD dimension: 6.7x6.9cm
Unit dimension: 13.1x10.1x4.4cm

Unit dimension: 11x10x5.7cm
LCD Dimension: 5.2cmx5.2cm

LCD dimension: 7.6x7.9cm
Unit dimension: 12.2x9.4x7.2cm

BLOOD PRESSURE MONITORS

ECONOMIC MODELS



Monitoring
Solutions



MODEL: AME DBP 1312

MODEL: AME 702V

TECHNOLOGY



FEATURES

Average last 3 results
Automatic power-off
Deluxe carry case
Backlight optional

2×60 memories with date and time
Talking optional
Bluetooth optional

2×99 memories with date and time

DESIGN ELEMENTS

Large display
4×AAA batteries or Type-C

5 Colors available: AC Adaptor

White finish
(colors available with MOQ)

SPECIFICATIONS

LCD dimension: 6.5x4.5cm
Unit dimension: 12.2x9.4x7.2cm

LCD dimension: 4.6x6.2cm
Unit dimension: 13.3x10.3x6.4cm



Monitoring
Solutions

BLOOD PRESSURE MONITORS ECONOMIC MODELS



MODEL: AME BP 809



MODEL: AME LT P70

TECHNOLOGY



Smart
DETECT Smart
INDICATE



Smart
CUFF Smart
DETECT Smart
INDICATE

FEATURES

Talking optional
Backlight Optional
Average last 3 results
Automatic Power Off

2×60 memories with date and time
Bluetooth optional

2×90 memories with date and time
Low battery detection
Optional Screen Color Changes with
BP Result: ■ ■ ■

DESIGN ELEMENTS

2 Colors available:
Large LED Display

4×AAA batteries or Type C or
AC adapter

4×AAA batteries or Type-C

SPECIFICATIONS

LCD dimension: 48×35mm
Unit dimension: 115×93×64mm

LCD dimension: 6.6×7cm
Unit dimension: 8.5×10.5×4.5cm

BLOOD GLUCOSE METERS

Welcome to a new era of diabetes management with AME's Blood Glucose Meters. Precision meets simplicity in our advanced devices, offering seamless blood glucose monitoring for enhanced well-being. Designed with user-friendly interfaces, AME's Glucose Meters ensure accurate readings at your fingertips, empowering individuals to take control of their health. AME's comprehensive solution integrates into your daily routine, making glucose monitoring a stress-free experience. Experience the future of diabetes



FEATURES

Our products are designed with a focus on accuracy, safety, and user convenience, ensuring reliable performance. Each device is manufactured to meet international quality and regulatory standards, providing healthcare professionals and consumers with confidence in every use.



DAILY ALARM

The daily alarm feature ensures timely reminders



MEASURING TIME 5 SECONDS

A very quick and accurate measuring time of just 5 seconds



STRIP EJECTOR

Allows for hygienic and hands-free disposal of used test strips



BLUETOOTH COMPATIBLE

Assay method used is compatible with Bluetooth. (works with GlucoJoy app)



BACK-LIGHT

Enhances visibility, enabling clear and easy use of the digital product even in low-light conditions



360 READING MEMORIES

The device has a reading memory of upto 360 readings



HYPO-HYPER ALARMS

The device beeps if the readings are outside of normal range



Monitoring
Solutions

AME BLOOD GLUCOSE METERS



AME BGM713r

AME BGM709b

TECHNOLOGY



Daily
ALARM



5
SECONDS



Strip
EJECTOR



Back-
LIGHT



360
Reading
MEMORIES



HYPO/HYPER
ALARMS



BLUETOOTH
Compatible

SPECIFICATIONS

Measurement range: 0.5–33.3 mmol/L (9–600 mg/dL)
 Sample: Fresh capillary whole blood, venous blood, neonatal blood
 Sample size: Approximate 0.6 microlitre
 Calibration: Plasma equivalent glucose values
 Battery life: Approximately 1,000 tests
 Unit of measure: mmol/L or mg/dL; switchable
 Memory: 360 blood glucose results with date and time
 Operating environment: 5°C45°C (41°F113°F); 10~90% RH (non-condensing)
 Meter storage environment: -20°C55°C (-4°F131°F); 10~95% RH (non-condensing)
 Strip storage environment: 1°C30°C (33.8°F86°F)
 Altitude: Up to 10,000 feet (3,048 meters) above sea level
 Hematocrit: 0%~70%
 Shelf life: Three years

Assay method: Glucose dehydrogenase biosensor

Battery: Two 3.0V CR2032 button battery, replaceable
 Size: 78*63.4*14mm(L*W*H)
 LCD size: 32*36.6mm(L*H)
 Weight: Approximate 40g, battery not included

Assay method: Glucose dehydrogenase biosensor with Bluetooth. (works with GlucoJoy app)
 Battery: 2×AAA alkaline batteries
 Size: 90*61*23.6 mm (L*W*H)
 LCD size: 35.4*30.9 mm (L*H)
 Weight: Approximate 51.6 g, battery not included

DIGITAL THERMOMETERS

Thermometers have come a long way since the Mercury filled counterparts . AME Thermometers are equipped with a large selection of cutting edge technologies which are reliable and effective in a home or clinical setting. Take look at the various models suitable for all applications.

CATEGORIES

-  Rigid Tip
-  Flexible Tip
-  Ear & Non-Contact



SMART TECHNOLOGIES

SMART TEMP

Pre-heated probe technology increases accuracy and reliability of temperature readings. It increases accuracy and stability by reducing the effect of temperature variations between the body and measurement site

SMART SELECT

Help customers select user age prior to measurement. Clinical Studies shows that children's body temperature is different from that of adults , hence AME SmartSelect Technology helps provide accurate measurement

SMART PREDICT

Employs the use of our patented algorithm to calculate the actual temperature in short time. upto 5S to 20S predictive reading

SMART POSITION

Guarantees temperature reliability through probe placement stability during the measurement process. The indicator light will pulse while the measurement is in progress. A steady light at the completion shows the probe has not moved

SMART PRECISION

The thermometer reading shows two decimals ensuring the necessary precision needed in clinical settings.

FEATURES



RIGID Tip



FLEXIBLE Tip



EAR Therm.



Non-CONTACT



1 SEC Reading



Switchable C INTO F



BLUE-Tooth



Automatic POWER OFF



Fever ALARM



TALKING



Distance Indicated LIGHT



Jumbo LCD



LED/LCD Display



Back-LIGHT



Monitoring
Solutions

DIGITAL THERMOMETERS RIGID TIP



MODEL: AME DMT 4161

FEATURES

- Rigid tip
- Fever alarm
- Last reading recall
- Dual scale with °C/°F
- Predictive measuring optional
- 10s/20s/30s response time
- Automatic power-off

DESIGN ELEMENTS

- Waterproof
- Beeps
- Low battery indicator
- White Finish

SPECIFICATIONS

- Unit Dimension: 14.3 x 2.5 x 1.4cm
- Unit Weight: 16grams (Battery Included)

FEATURES



**RIGID
TIP**



**SMART
PREDICT**



**SWITCHABLE
C INTO F**



**BACK-
LIGHT**



**FEVER
ALARM**



**AUTOMATIC
POWER OFF**

DIGITAL THERMOMETERS

RIGID TIP



Monitoring Solutions



MODEL: AME DMT 4759

MODEL: AME DMT 4763

TECHNOLOGY



Jumbo LCD



Rigid TIP



Switchable C INTO F



Back-LIGHT



Fever ALARM



Automatic POWER OFF



Smart PREDICT

FEATURES

- Rigid tip
- Fever alarm
- Last reading recall
- Dual scale with °C/°F
- 10s/20s/30s response time
- Automatic power-off

10 Reading Memories

Predictive measuring optional

DESIGN ELEMENTS

- Waterproof
- Beeps
- Low battery indicator
- White Finish

Jumbo LCD Size



Monitoring Solutions

DIGITAL THERMOMETERS

RIGID TIP



MODEL: AME DMT 4119

MODEL: AME DMT 101

MODEL: AME DMT 412

TECHNOLOGY



FEATURES

Rigid tip
Fever alarm
Last reading recall
Dual scale with °C/°F
Automatic power-off

10s/20s/30s response time
Fever line indicator (Optional)

10s/20s/30s response time

DESIGN ELEMENTS

Beeps
White Finish

Jumbo LCD

DIGITAL THERMOMETERS

FLEXIBLE TIP



Monitoring Solutions



MODEL: AME DMT 4336

MODEL: AME DMT 4760

MODEL: AME DMT 4765

TECHNOLOGY

 Flexible TIP	 Switchable C INTO F	 Flexible TIP	 Switchable C INTO F	 Blue-TOOTH	 Smart PREDICT	 Flexible TIP	 Switchable C INTO F	 Smart PREDICT
 Automatic POWER OFF	 Fever ALARM	 Automatic POWER OFF	 Fever ALARM	 Jumbo LCD	 Smart PRECISION	 Automatic POWER OFF	 Fever ALARM	 Back-LIGHT

FEATURES

Beeps
Flexible tip
Fever alarm
Waterproof
Last reading recall
Dual scale with °C/°F
10s/20s/30s response time
Automatic power-off

Beeps
Flexible Tip
Fever Alarm
Waterproof
Bluetooth Optional
10 reading memory
Predictive Measuring Optional
Jumbo LCD
High Precision
10s/20s/30c Options
Dual Scale C & F Convertible
Low Battery Indicator
Automatic Power Off

Beeps
Flexible Tip
Fever Alarm
Waterproof
Last Reading Recall
Backlight Optional
10 Reading Memories
Predictive Measuring Optional
10s/20s/30c Options
Dual Scale C & F Convertible
Low Battery Indicator
Automatic Power Off



Monitoring
Solutions

DIGITAL THERMOMETERS FLEXIBLE TIP



MODEL: AVE DMT 437 Cartoon Series

TECHNOLOGY



Flexible
TIP



Switchable
C INTO F



Automatic
POWER OFF



Fever
ALARM

FEATURES

- Flexible Tip
- Waterproof
- Fever alarm
- Last reading recall
- Dual scale with °C/°F
- 10s/20s/30s response time
- Automatic power-off

DESIGN ELEMENTS

- Beeps
- Colorful designs

DIGITAL THERMOMETERS

EAR & NON-CONTACT



Monitoring
Solutions



MODEL:
AME DET 1215



MODEL:
AME DET 1115

TECHNOLOGY



SMART TEMP

Pre-heated probe technology increases accuracy and reliability of temperature readings. It increases accuracy and stability by reducing the effect of temperature variations between the body and measurement site



SMART POSITION

Guarantees temperature reliability through probe placement stability during the measurement process. The indicator light will pulse while the measurement is in progress. A steady light at the completion shows the probe has not moved

FEATURES

- Measures in ear
- Pre-heated technology
- Preheating time 6-10s
- 30 Memories with date & time
- 3 seconds reading
- Dual scale with °C/°F
- Automatic power-off
- Bluetooth optional
- 60s timer optional
- With probe covers

DESIGN ELEMENTS

- Beeps; facilitates precise timing with audible cues every 15 seconds to 60 seconds
- Indicator light
- Black Finish
- Rechargeable battery optional
- Backlight optional
- Battery for AME DET 1215: 1900mAh li-battery
- Battery for AME DET 1115: 2*AA battery

STANDARD AME ACCESSORIES



Probe Cover



High-end recharge-
able storage box



Monitoring
Solutions

DIGITAL THERMOMETERS EAR & NON-CONTACT



MODEL: AME F7

MODEL: AME DET 3023

TECHNOLOGY



FEATURES

Dual scale with °C/°F
1 second reading

Fever alarm
9 group memories
Mute switching function
Disposable probe cover to avoid cross infection
Ejection of probe covers
Cover detection function to avoid misoperation

Measurement of forehead, objects & Environment
30 reading memories
Talking optional
Bluetooth optional
Automatic power-off

DESIGN ELEMENTS

Beeps

White finish
Large LED display
With probe covers

Black finish
Replaceable battery

DIGITAL THERMOMETERS

EAR & NON-CONTACT



Monitoring
Solutions



TECHNOLOGY



SMART POSITION

Pre-heated probe technology increases accuracy and reliability of temperature readings. It increases accuracy and stability by reducing the effect of temperature variations between the body and measurement site

MODEL: AME AET U613A

FEATURES

- Measurement of forehead, objects & Environment
- No probe covers needed
- 10 reading memories
- 1 second reading
- Bluetooth optional
- Dual scale with °C/°F
- Automatic power-off


DESIGN ELEMENTS

- Compact ergonomic size
- Beeps
- Backlight optional
- Talking optional
- Replaceable battery
- Colors available



CLINICAL ESSENTIALS

Elevate your clinical practice with AME's precision-engineered Stethoscopes and Sphygmomanometers. Crafted with a blend of advanced technology and ergonomic design, our stethoscopes deliver clear, accurate auscultation for healthcare professionals. Paired with our Sphygmomanometers, featuring precision measurement and durability, AME ensures reliable healthcare monitoring. Trust in the superior quality and performance of AME's medical devices, designed to meet the exacting standards of healthcare providers worldwide.



PRODUCT LINE

Our products are designed with a focus on accuracy, safety, and user convenience, ensuring reliable performance. Each device is manufactured to meet international quality and regulatory standards, providing healthcare professionals and consumers with confidence in every use.



CLASSIC STETHOSCOPE

Dual head Professional stethoscope with rugged durability and high acoustic response



ANEROID SPHYGMOMANOMETER

Aneroid Sphygmomanometer with Free Stethoscope and Free Calibration Key



DIPLOMAT STETHOSCOPE

Dual head Nurses stethoscope, suitable for home use




DIGITAL SPHYGMOMANOMETER

Digital palm type, mercury free sphygmomanometer

CLINICAL ESSENTIALS STETHOSCOPES



Monitoring
Solutions



MODEL: AME 3006

MODEL: AME 3002



FEATURES

Dual head type
Rugged Durability
High Acoustic Response

Chest piece is machined from surgical stainless steel

Aluminum Stethoscope

DESIGN ELEMENTS

Overall length 770 mm

Foam lined box
Two extra sets of soft sealing eartips
One extra diaphragm



Monitoring
Solutions

CLINICAL ESSENTIALS SPHYGMOMANOMETERS



AME 2001-3001

AME DELUX PLUS



AME 1019

**AME DIGITAL PALM TYPE
MERCURY FREE SPHYGMOMANOMETER**



**ANEROID
Sphygmomanometer**



**DIGITAL
Sphygmomanometer**

FEATURES

300mmHg, no pin stop manometer
Available in five sizes
Comes with standard air release valve
Free Calibration Key Included
Free Stethoscope included

Measuring range: Pressure: (0-40)kPa, (0-300)mmHg
Comes with Armband & Digital Inflation bulb
Expected service life of sphygmomanometer is 3years

DESIGN ELEMENTS

Black color
Material is Nylon
Packed into black vinyl zipper case

Resolution: Liquid crystal column: 0.26kPa(2mmHg)
Digital display: 0.13kPa(1mmHg)
Power supply: (AAA) * 2 alkaline dry battery
Suitable Arm circumference: 220mm-320mm
Displays Systolic, Diastolic and Pulse with special technique.

WEIGHING SCALES

Discover the precision of AME's Weighing Scales, seamlessly blending everyday functionality with cutting-edge technology. Our Scales offer reliability and accuracy, ensuring a fundamental tool for health-conscious individuals. Our latest models are equipped with Bluetooth connectivity for a comprehensive health summary.

CATEGORIES

-  Smart Scales
-  Digital Glass Scales



TECHNOLOGIES

Our products are designed with a focus on accuracy, safety, and user convenience, ensuring reliable performance. Each device is manufactured to meet international quality and regulatory standards.

LED DISPLAY

Bright, high-visibility numbers for quick viewing in any light

BLUETOOTH COMPATIBLE

Syncs seamlessly with your smartphone app for tracking and analysis

TOUGHENED GLASS

Durable, shatter-resistant platform for long-lasting use

LCD DISPLAY

Clear, easy-to-read screen for precise weight readings

KG/LB CONVERTIBLE

Switch instantly between kilograms and pounds for flexible measurement

AUTOMATIC ON / OFF

Powers on when stepped on and off when not in use to save energy

SMART SCALE FEATURES

 WEIGHT	 FAT RATIO	 BMI	 VISCERAL FAT	 BMR	 FAT CONTROL	 FAT MASS	 FAT-FREE MASS
 BONE MASS	 MUSCLE MASS	 MUSCLE RATE	 MUSCLE CONTROL	 SKELETAL MUSCLE RATE	 SKELETAL MUSCLE MASS	 STANDARD WEIGHT	 WEIGHT CONTROL
 BODY TYPE	 OBESITY LEVEL	 HEALTH EVALUATION	 BODY AGE	 HEALTH SCORE	 PROTEIN RATIO	 WATER RATIO	 SUBCUTANEOUS FAT RATIO

WEIGHING SCALES

SMART SCALES



Monitoring Solutions



AME CF577BLE+WIFI 8 Electrode Dual Frequency Smart Body Fat Scale

FEATURES

- User Selection
- Works with WIFI
- Works with Bluetooth
- Easy to use AME App
- BIA Measurement
- Kg/Lb Convertible
- 180kg Capacity

DESIGN ELEMENTS

- 3.5 Inch TFT color large display
- 81.5mm*53.5mm Size
- Sturdy Platform
- 5mm Thickness with toughened glass
- Power: 3*AAA Batteries

SPECIFICATIONS

Scale Size: 335*335*35mm



TFT Color Display



User Selection/Profile



Guided Use

FEATURES

- WEIGHT
- BONE MASS
- BODY TYPE
- FAT RATIO
- MUSCLE MASS
- OBESITY LEVEL
- BMI
- MUSCLE RATE
- HEALTH EVALUATION
- VISCERAL FAT
- MUSCLE CONTROL
- BODY AGE
- BMR
- SKELETAL MUSCLE RATE
- HEALTH SCORE
- FAT CONTROL
- SKELETAL MUSCLE MASS
- PROTEIN RATIO
- FAT MASS
- STANDARD WEIGHT
- WATER RATIO
- FAT-FREE MASS
- WEIGHT CONTROL
- SUBCUTANEOUS FAT RATIO
- BLUETOOTH Compatible



Monitoring Solutions

WEIGHING SCALES SMART SCALES



AME CF539 Digital Body Fat Scale with Solar Charging

FEATURES

- Solar/Light Charging
- Kg/Lb Convertible
- 180kg Capacity
- Easy to use AME App
- App provides all the listed icon features

DESIGN ELEMENTS

- Spectacular White Design
- Sturdy Platform
- 5mm Thickness with toughened glass
- Power: 3*AAA Batteries

SPECIFICATIONS

Scale Size: 280*280*23mm

FEATURES



WEIGHT



BONE MASS



BODY TYPE



FAT RATIO



MUSCLE MASS



OBESITY LEVEL



BMI



MUSCLE RATE



HEALTH EVALUATION



VISCERAL FAT



MUSCLE CONTROL



BODY AGE



BMR



SKELETAL MUSCLE RATE



HEALTH SCORE



FAT CONTROL



SKELETAL MUSCLE MASS



PROTEIN RATIO



FAT MASS



STANDARD WEIGHT



WATER RATIO



FAT-FREE MASS



WEIGHT CONTROL



SUBCUTANEOUS FAT RATIO



BLUETOOTH Compatible

WEIGHING SCALES

DIGITAL GLASS SCALES



Monitoring Solutions



AME CW 278



AME CW 286



AME CW 310

TECHNOLOGY



FEATURES

180 kg Weighing Capacity
 Kg/Lb Convertible
 Sturdy Platform
 5mm Thickness with toughened glass
 Power: 3*AAA Batteries

LCD Display
 Apple White Finish
 Blue Backlight Display
 Scale Size: 280*280*24mm

Hidden LED Display
 Metallic Powder Finish
 Scale Size: 280*280*24mm

Hidden LED Display
 Optional Designs
 Scale Size: 260*260*24mm



Monitoring
Solutions

WEIGHING SCALES DIGITAL GLASS SCALES



AME CW 516BLE

AME CF 598 BL

TECHNOLOGY



FEATURES

LCD Display
Variable Colors and Designs
180 kg Weighing Capacity
Kg/Lb Convertible
Blue Backlight Display
Scale Size: 280*280*24mm
Sturdy Platform
5mm Thickness with toughened glass
Power: 3*AAA Batteries
Bluetooth Version: BLE4.0, Works with iOS8.0&Android4.3 OR Higher Version
Measurements: 24 measures, include Body weight, Fat, Water, BMI, Bone, Muscle, Visceral Fat, BMR, etc
Sensors: 4x Sensors - High Sensitive
AUTO Function: Step on, Auto-off;
Smart Reminder: Low Battery Reminder, Overload Reminder.
Free APP "Unique Health" or Customized APP or provide the SDK for customization

Respiratory Care Solutions



Breathing comfort and respiratory health are at the heart of our care philosophy. AME's Respiratory Care Solutions provide reliable, clinically tested devices that help relieve congestion, support therapy for respiratory conditions, and improve overall air quality.

Our range includes nebulizers, steam inhalers, and humidifiers, all designed for ease of use, quiet operation, and consistent performance. Whether for professional clinical settings or home care, these solutions deliver dependable support for healthier, easier breathing every day.

Product Range



Nebulizers



Steam Inhalers



Humidifiers



NEBULIZERS

Nebulizers are Class IIb Medical Devices which are suitable for the treatment of Asthma and their symptoms. AME Nebulizers utilize a wide variety of innovative technologies which are beneficial in home settings or clinical settings. Take a look at our models suitable for all applications.

CATEGORIES

-  Compressor Nebulizer - **CLASSIC**
-  Compressor Nebulizer - **CLASSIC BOLD**
-  Compressor Nebulizer - **ECO SERIES**
-  Compressor Nebulizer - **ANIMAL SERIES**
-  Mesh Nebulizer - **LATEST**



TECHNOLOGIES

Our products are designed with a focus on accuracy, safety, and user convenience, ensuring reliable performance. Each device is manufactured to meet international quality and regulatory standards, providing healthcare professionals and consumers with confidence in every use.



DUAL SPEED
2 speed chamber allows for adjustment of mist produced to help save medicines & avoid choking hazard for infants



NON-STOP USE
Continuous use model which can be used day & night for more than 100 hours without stopping; suitable for home use small clinics



LESS NOISE
Strong motor for long running cycles with less noise, providing comfortable environment



CYCLIC USE
Suitable for smaller families when used intermittently - 40 minutes continuous use followed by 10 mins rest



USB-C TYPE
To provide adaptability across borders, we utilize the USB C to ensure continuity in therapy when travelling



WARRANTY 5 YEARS
5 year warranty on manufacturing defects



NOISE LESS
Perfect for situations where noise will disturb therapy



RECHARGEABLE TYPE
No more running out of batteries as the product is re-chargeable.

COMPRESSOR NEBULIZER CLASSIC



Respiratory
Solutions



FEATURES



MODEL: AME MCN S600B

SPECIFICATIONS

Power:	AC 230V / 50Hz
Power Consumption:	180VA
Nebulizer medication Capacity:	8ml
Partical Size:	Approx.(75.74±0.02)%; within 0.5 to 5um range
MMAD:	2.81±0.14µm
Average Nebulization Rate:	≥0.2ml/min
Noise Level:	Below 65dBA (DIN EN 13544-1)
Compressor Pressure Range:	30 ~ 50 Psi / 206Kpa~345Kpa
Operation Pressure Range:	9.5 ~ 20 Psi / 65.5Kpa~137Kpa
Liter Flow Range:	6~10 lpm
Operation Temperature Range:	10°C ~40°C (50°F ~104°F)
Operation Humidity Range:	10% ~95% RH
Storage Temperature Range:	-20°C ~70°C (-4°F~158°F)
Storage Humidity Range:	10%~95% RH
Atmospheric pressure	500~1060hPa

STANDARD AME ACCESSORIES



2 Speed
Nebulizer Chamber



Adult Mask



Child Mask



Air Tube



Mouth
Piece



NosePiece
(Optional)



Respiratory
Solutions

COMPRESSOR NEBULIZER CLASSIC



AME MCNS600A



AME MCN S600E

TECHNOLOGY



DUAL
Speed



NON-STOP
Use



LESS
Noise



5 YEAR
WARRANTY

SPECIFICATIONS

Power:	AC 230V / 50Hz
Power Consumption:	180VA
Nebulizer medication Capacity:	8ml
Partical Size:	Approx.(75.74±0.02)%; within 0.5 to 5um range
MMAD:	2.81±0.14µm
Average Nebulization Rate:	≥0.2ml/min
Noise Level:	Below 65dBA (DIN EN 13544-1)
Compressor Pressure Range:	30 ~ 50 Psi / 206Kpa~345Kpa
Operation Pressure Range:	9.5 ~ 20 Psi / 65.5Kpa~137Kpa
Liter Flow Range:	6~10 lpm
Operation Temperature Range:	10°C ~40°C (50°F ~104°F)
Operation Humidity Range:	10% ~95% RH
Storage Temperature Range:	-20°C ~70°C (-4°F~158°F)
Storage Humidity Range:	10%~95% RH
Atmospheric pressure	500~1060hPa

STANDARD AME ACCESSORIES



2 Speed
Nebulizer Chamber



Adult Mask



Child Mask



Air Tube



Mouth
Piece



NosePiece
(Optional)

COMPRESSOR NEBULIZER

CLASSIC BOLD



Respiratory
Solutions



FEATURES



MODEL: AME MCN S600G

SPECIFICATIONS

Power:	AC 230V / 50Hz
Power Consumption:	180VA
Nebulizer medication Capacity:	8ml
Partical Size:	Approx.(75.74±0.02)% within 0.5 to 5um range
MMAD:	2.81±0.14µm
Average Nebulization Rate:	≥0.2ml/min
Noise Level:	Below 65dBA (DIN EN 13544-1)
Compressor Pressure Range:	30 ~ 50 Psi / 206Kpa~345Kpa
Operation Pressure Range:	9.5 ~ 20 Psi / 65.5Kpa~137Kpa
Liter Flow Range:	6~10 lpm
Operation Temperature Range:	10°C ~40°C (50°F ~104°F)
Operation Humidity Range:	10% ~95% RH
Storage Temperature Range:	-20°C ~ 70°C (-4°F~158°F)
Storage Humidity Range:	10%~95% RH
Atmospheric pressure	500~1060hPa

STANDARD AME ACCESSORIES





Respiratory Solutions

COMPRESSOR NEBULIZER CLASSIC BOLD



AME MCNS600 F



AME MCN S600 MD

TECHNOLOGY



DUAL Speed



NON-STOP Use



LESS Noise



5 YEAR WARRANTY

SPECIFICATIONS

Power:	AC 230V / 50Hz
Power Consumption:	180VA
Nebulizer medication Capacity:	8ml
Partical Size:	Approx.(75.74±0.02)% within 0.5 to 5um range
MMAD:	2.81±0.14µm
Average Nebulization Rate:	≥0.2ml/min
Noise Level:	Below 65dBA (DIN EN 13544-1)
Compressor Pressure Range:	30 ~ 50 Psi / 206Kpa~345Kpa
Operation Pressure Range:	9.5 ~ 20 Psi / 65.5Kpa~137Kpa
Liter Flow Range:	6~10 lpm
Operation Temperature Range:	10°C ~40°C (50°F ~104°F)
Operation Humidity Range:	10% ~95% RH
Storage Temperature Range:	-20°C ~ 70°C (-4°F~158°F)
Storage Humidity Range:	10%~95% RH
Atmospheric pressure	500~1060hPa

STANDARD AME ACCESSORIES



2 Speed Nebulizer Chamber



Adult Mask



Child Mask



Air Tube



Mouth Piece



Nose Piece (Optional)

COMPRESSOR NEBULIZER

ECO SERIES



Respiratory
Solutions



FEATURES



MODEL: AME MCN S600MN

SPECIFICATIONS

Power:	AC 120V 60Hz, 1.3A or 230V 50Hz 0.8A
Power Consumption:	230V 50Hz 0.8A
Nebulizer medication Capacity:	6 ml
Partical Size:	≥65% within 0.5 ~5um range
Average Nebulization Rate:	≥0.2 ml/min
Noise Level:	≤65dBA(acc.to DIN EN 13544-1 section 26)
Compressor Pressure Range:	30 ~ 50 Psi / 206Kpa~345Kpa
Operation Pressure Range:	Approx.7.0~13.0 Psi/48.0~89.0Kpa
Liter Flow Range:	Approx.3.0~6.0 lpm
Operation Temperature Range:	10 ~40°C, 30 ~ 85 %
Operation Humidity Range:	10% ~95% RH
Storage Temperature Range:	-20°C ~ 70°C (-4°F~158°F)
Storage Humidity Range:	20 ~ 75 % RH
Atmospheric pressure	800~1060hPa
Expected service life	Approx.1000H
Working and rest time	30 mins On, 30 minsOFF
Safety level	Type BF ClassII

STANDARD AME ACCESSORIES





Respiratory
Solutions

COMPRESSOR NEBULIZER ECO SERIES



AME MCN S600MA



AME MCN S600ME

TECHNOLOGY



DUAL
Speed



NON-STOP
Use



CYCLIC
Use



5 YEAR
WARRANTY

SPECIFICATIONS

Power:	AC 230V / 50Hz
Power Consumption:	180VA
Nebulizer medication Capacity:	8ml
Partical Size:	Approx.(75.74±0.02)% within 0.5 to 5um range
MMAD:	2.81±0.14µm
Average Nebulization Rate:	≥0.2ml/min
Noise Level:	Below 65dBA (DIN EN 13544-1)
Compressor Pressure Range:	30 ~ 50 Psi / 206Kpa~345Kpa
Operation Pressure Range:	9.5 ~ 20 Psi / 65.5Kpa~137Kpa
Liter Flow Range:	6~10 lpm
Operation Temperature Range:	10°C ~40°C (50°F ~104°F)
Operation Humidity Range:	10% ~95% RH
Storage Temperature Range:	-20°C ~70°C (-4°F~158°F)
Storage Humidity Range:	10%~95% RH
Atmospheric pressure	500~1060hPa

STANDARD AME ACCESSORIES



2 Speed
Nebulizer Chamber



Adult Mask



Child Mask



Air Tube



Mouth
Piece



NosePiece
(Optional)

COMPRESSOR NEBULIZER

ECO SERIES



Respiratory
Solutions



AME MCN S600MC



AME AN106

TECHNOLOGY



DUAL
Speed



NON-STOP
Use



CYCLIC
Use



5 YEAR
WARRANTY

SPECIFICATIONS

Power:	AC 230V / 50Hz
Power Consumption:	180VA
Nebulizer medication Capacity:	8ml
Partical Size:	Approx.(75.74±0.02)% within 0.5 to 5um range
MMAD:	2.81±0.14µm
Average Nebulization Rate:	≥0.2ml/min
Noise Level:	Below 65dBA (DIN EN 13544-1)
Compressor Pressure Range:	30 ~ 50 Psi / 206Kpa~345Kpa
Operation Pressure Range:	9.5 ~ 20 Psi / 65.5Kpa~137Kpa
Liter Flow Range:	6~10 lpm
Operation Temperature Range:	10°C ~40°C (50°F ~104°F)
Operation Humidity Range:	10% ~95% RH
Storage Temperature Range:	-20°C ~70°C (-4°F~158°F)
Storage Humidity Range:	10%~95% RH
Atmospheric pressure	500~1060hPa

STANDARD AME ACCESSORIES





Respiratory
Solutions

COMPRESSOR NEBULIZER ECO SERIES



AME JLN 361 AIR



AME JLN 2300 AIR

TECHNOLOGY



SPECIFICATIONS

Power:	AC 230V / 50Hz
Power Consumption:	180VA
Nebulizer medication Capacity:	8ml
Partical Size:	Approx.(75.74±0.02)% within 0.5 to 5um range
MMAD:	2.81±0.14µm
Average Nebulization Rate:	≥0.2ml/min
Noise Level:	Below 60dBA (DIN EN 13544-1)
Compressor Pressure Range:	30 ~ 50 Psi / 206Kpa~345Kpa
Operation Pressure Range:	9.5 ~ 20 Psi / 65.5Kpa~137Kpa
Liter Flow Range:	6~10 lpm
Operation Temperature Range:	10°C ~40°C (50°F ~104°F)
Operation Humidity Range:	10% ~95% RH
Storage Temperature Range:	-20°C ~70°C (-4°F~158°F)
Storage Humidity Range:	10%~95% RH
Atmospheric pressure	500~1060hPa

STANDARD AME ACCESSORIES



2 Speed
Nebulizer Chamber



Adult Mask



Child Mask



Air Tube



Mouth
Piece



Nose Piece
(Optional)

COMPRESSOR NEBULIZER

ECO SERIES



Respiratory
Solutions



FEATURES



MODEL: AME 303

SPECIFICATIONS

Power:	AC 120V 60Hz, 1.3A or 230V 50Hz 0.8A
Power Consumption:	230V 50Hz 0.8A
Nebulizer medication Capacity:	6 ml
Partical Size:	≥65% within 0.5 ~5um range
Average Nebulization Rate:	≥0.2 ml/min
Noise Level:	≤65dBA(acc.to DIN EN 13544-1 section 26)
Compressor Pressure Range:	30 ~ 50 Psi / 206Kpa~345Kpa
Operation Pressure Range:	Approx.7.0~13.0 Psi/48.0~89.0Kpa
Liter Flow Range:	Approx.3.0~6.0 lpm
Operation Temperature Range:	10 ~40°C, 30 ~ 85 %
Operation Humidity Range:	10% ~95% RH
Storage Temperature Range:	-20°C ~ 70°C (-4°F~158°F)
Storage Humidity Range:	20 ~ 75 % RH
Atmospheric pressure	800~1060hPa
Expected service life	Approx.1000H
Working and rest time	30 mins On, 30 minsOFF
Safety level	Type BF ClassII
Colors available	White, Blue

STANDARD AME ACCESSORIES





Respiratory
Solutions

COMPRESSOR NEBULIZER

DC TYPE SERIES



FEATURES



MODEL: AME NB505

SPECIFICATIONS

Charging Time	60 minutes
Battery capacity	1200 mAh
Charging attribute	12v 1.0A Max
Residual amount	≤0.8mL
Output gas flow	>5ML/Min
Median particle size	3.5um±25%
Working noise	≤50dB(A)
Power	12V Power Adapter Included
Functionality	One Touch key Function
Digital Display	5.3 Inch Screen

STANDARD AME ACCESSORIES



COMPRESSOR NEBULIZER

ANIMAL SERIES



Respiratory
Solutions



FEATURES



MODEL: AMEMCNS600 J

SPECIFICATIONS

Power:	AC 230V / 50Hz
Power Consumption:	180VA
Nebulizer medication Capacity:	8ml
Partical Size:	Approx.(75.74±0.02)% within 0.5 to 5um range
MMAD:	2.81±0.14µm
Average Nebulization Rate:	≥0.2ml/min
Noise Level:	Below 65dBA (DIN EN 13544-1)
Compressor Pressure Range:	30 ~ 50 Psi / 206Kpa~345Kpa
Operation Pressure Range:	9.5 ~ 20 Psi / 65.5Kpa~137Kpa
Liter Flow Range:	6~10 lpm
Operation Temperature Range:	10°C ~40°C (50°F ~104°F)
Operation Humidity Range:	10% ~95% RH
Storage Temperature Range:	-20°C ~70°C (-4°F~158°F)
Storage Humidity Range:	10%~95% RH
Atmospheric pressure	500~1060hPa

STANDARD AME ACCESSORIES



2 Speed
Nebulizer Chamber



Adult Mask



Child Mask



Air Tube



Mouth
Piece



NosePiece
(Optional)



Respiratory
Solutions

COMPRESSOR NEBULIZER ANIMAL SERIES



AME MCN S600 H



AME MCN S600 MB

TECHNOLOGY



DUAL
Speed



NON-STOP
Use



LESS
Noise



5 YEAR
WARRANTY

SPECIFICATIONS

Power:	AC 230V / 50Hz
Power Consumption:	180VA
Nebulizer medication Capacity:	8ml
Partical Size:	Approx.(75.74±0.02)% within 0.5 to 5um range
MMAD:	2.81±0.14µm
Average Nebulization Rate:	≥0.2ml/min
Noise Level:	Below 65dBA (DIN EN 13544-1)
Compressor Pressure Range:	30 ~ 50 Psi / 206Kpa~345Kpa
Operation Pressure Range:	9.5 ~ 20 Psi / 65.5Kpa~137Kpa
Liter Flow Range:	6~10 lpm
Operation Temperature Range:	10°C ~40°C (50°F ~104°F)
Operation Humidity Range:	10% ~95% RH
Storage Temperature Range:	-20°C ~70°C (-4°F~158°F)
Storage Humidity Range:	10%~95% RH
Atmospheric pressure	500~1060hPa

STANDARD AME ACCESSORIES



2 Speed
Nebulizer Chamber



Adult Mask



Child Mask



Air Tube



Mouth
Piece



Nose Piece
(Optional)

MESH NEBULIZER

LATEST



Respiratory
Solutions



FEATURES



CYCLIC
USE



USB-C
TYPE



RECHARGEABLE
TYPE



NOISE
LESS

MODEL: AME MESH600 MF

SPECIFICATIONS

Method of Operation	Active Vibrating Mesh Technology
Power Supply	10mins Auto off, Internal power supply, Lithium battery only
Power Consumption	Approx. 2.0W
Auto Off	10 mins Auto off (with rechargeable battery)
Light Indications	Blue Light: During operation
	Orange Light: Battery is Low power
	Orange Light Blinking: Battery is Charging
	Blue Light Blinking: Battery has been fully Charged
Vibrating Frequency	Approx. 107 kHz +/- 10%
Nebulization Rate	≥ 0.25 ml/min
Particle Size	Approx. (74.53 \pm 3.01)% within 0.5 to 5 μ m range
MMAD	2.93 \pm 0.19 μ m
Capacity of Medication Cup	Max 6ml or 12ml
Noise Level	≤ 50 dBA (1 meter distance)
Battery Charging Voltage and Current	5.0V, Max. 2A
Battery Charging Time	Approx. 180 mins
Battery Non-stop Usage Time	Approx. 120 mins
Lithium Titanate Battery Cycle Life	The cycle time is not less than 300 times

STANDARD AME ACCESSORIES



Adult Mask



Child Mask



USB C Cable



Respiratory
Solutions

MESH NEBULIZER LATEST



AME MESH600 G



AME MESH600 C



AME MESH600 E

SPECIFICATIONS

Method Of Operation:	Active Vibrating Mesh Technology
Power Supply:	Lithium Titanate Battery (DC 4.8V Rechargeable)
Power Consumption:	Approx.2.0W
Auto Off:	10 Mins
Vibrating Frequency:	Approx.107 kHz +/- 10%
Nebulization Rate:	≥0.25 MI/Min
Particle Size:	Approx.(88.08+1.91)% within 0.5 to 5um range
MMAD:	2.25+0.35um
Capacity Of Medication Cup:	Max 12ml
Noise Level:	≤50dBA (1 Meter Distance)
Battery Chargeing Voltage And Current:	5.0V 1A~2A
Battery Chargeing Time:	Approx.90 Mins
Battery Non-Stop Usage Time:	Approx.100 Mins
Lithium Titanate Battery Cycle Life:	Cycle Life:0.5C Approx.7000 Times, 1C Approx.4000 Times
Weight:	Approx.110g
Operating Temperature And Humidity:	10 ~40°C, 30 ~ 85 % RH, 800~1060hPa
Storage And Delivery Condition: -	20~ 60°C, 20~ 75 % RH, 800~1060hPa
Safety Level:	Type B Class II

FEATURES



CYCLIC
Use



USB-C
TYPE



RECHARGE-
ABLE Type



NOISE
Less

STANDARD AME ACCESSORIES



Adult Mask



Child Mask



STEAM INHALERS

Introducing our latest innovation in respiratory care - the AME Steam Inhalers. Harnessing the power of steam therapy, it provides a soothing and effective solution for congestion, cough, and other respiratory issues. Elevate your well-being with the 2 types of AME Steam Inhalers, where advanced technology meets uncompromising quality for a breath of fresh, revitalizing air.

CATEGORIES

- AME EcoSteam
- AME HandySteam



FEATURES

Our products are designed with a focus on accuracy, safety, and user convenience, ensuring reliable performance. Each device is manufactured to meet international quality and regulatory standards, providing healthcare professionals and consumers with confidence in every use.



Maintains the temperature of the device when it exceeds limits, preventing damage.



Performs quick action by boiling water within 3 minutes



Works at 3 different settings 0/1/2 Control



Ensures easy cleanup after usage, also prevents particles from sticking to the surface



Allows easy storage & convenient portability, enhancing the usage



Respiratory Solutions

STEAM INHALERS Facial Steamer & Steam Inhaler



AME ECOSTEAM



AME HANDYSTEAM

TECHNOLOGY



FEATURES

PTC fast heater
Boiling within 3-min
Over-heating safety protection

Steam Control design
3 settings 0/1/2 Control

Detachable power cord
Prevent leakage design

DESIGN ELEMENTS

220-240VAC 50/60Hz 100W
LED power indicator

With facial mask, inhaler and cup
Capacity: 50ml

With facial mask and cup
ON/OFF Switch
Capacity: 30ml

HUMIDIFIER

Another one of our latest category in respiratory care is the AME Humidifiers. Crafted with customer usage and designed for optimal user experience, our humidifiers deliver targeted relief for respiratory discomfort. Harnessing the power of Hi-Tech Atomizers, it provides a soothing and effective solution for environments which have very dry air. Elevate your well-being with the AME Steam Humidifiers, where advanced technology meets uncompromising quality for a breath of fresh, revitalizing air. Breathe better, live healthier for a Better life.



FEATURES

Our products are designed with a focus on accuracy, safety, and user convenience, ensuring reliable performance. Each device is manufactured to meet international quality and regulatory standards, providing healthcare professionals and consumers with confidence in every use.



LARGE CAPACITY

Comes with a large 2 Litre capacity, ensuring to run for a long time



FAST ACTION

Performs quick action by ultra fast atomization of water



NIGHT LIGHT

Provides a soft rainbow color changing light, perfect at night



EASY TO CLEAN

The top lid comes off to allow easy wiping & cleaning



Respiratory Solutions

HUMIDIFIER Latest



FEATURES



LARGE CAPACITY
Comes with a large 2 Litre capacity, ensuring to run for a long time



FAST ACTION
Performs quick action by ultra fast atomization of water



EASY TO CLEAN
The top lid comes off to allow easy wiping & cleaning



NIGHT LIGHT
Provides a soft rainbow color changing light, perfect at night

MODEL: AME HMD01

FEATURES

- Ultra Fine and Fast Atomization
- Ambient Light Color Change
- Blue Light Control
- Non-stick pan technology
- Over-heating safety protection
- Suitable for all night usage.

DESIGN ELEMENTS

- 220-240VAC 50/60Hz 100W
- Capacity: 2000ml
- LED power indicator
- Blue Light Control
- 2 Colors available:



Rapid Diagnostic Solutions



Early detection leads to better decisions, and at AME we make it simple. Our Rapid Diagnostic Solutions are designed to deliver fast, accurate, and dependable results, empowering both healthcare professionals and individuals to make informed choices with confidence.

From pregnancy and fertility test kits to other in-vitro diagnostic tools, each product combines precision technology with ease of use, ensuring laboratory-quality accuracy in every setting. Reliable, certified, and convenient, AME brings trusted diagnostics to your fingertips.

Product Range



Pregnancy Test Kits



Fertility Test Kits



Drug of Abuse Kits

For More categories, download the complete IVD catalogue

PRODUCT LINE

- Pregnancy Test Kits
- Ovulation Tests
- Menopause Tests
- Drug of Abuse Test Kits
- Veterinary Test Kits
- Home Use Tests



PREGNANCY & FERTILITY TEST KITS

Welcome to the forefront of personalized diagnostics with AME's In-Vitro Diagnostics Kits. Our Pregnancy Test Kits redefine convenience and accuracy in home testing. Designed with precision, AME empowers users with reliable results, enabling proactive health management. Experience a new era in at-home healthcare with AME's In-Vitro Diagnostics Kits - where accuracy meets convenience.

CATEGORIES

- Mid Stream Type
- Cassette Type



FEATURES

Our products are designed with a focus on accuracy, safety, and user convenience, ensuring reliable performance. Each device is manufactured to meet international quality and regulatory standards, providing healthcare professionals and consumers with confidence in every use.



HOME USE
Provides convenience and accurate results in the comfort of your home



MIDSTREAM TYPE
Offers a convenient design, allowing users to directly test in the urine stream



FERTILITY DATE
Estimates the likely date of conception helping to understand and track pregnancy timeline



QUICK & ACCURATE
Ensures reliable results in minimal time, enhancing efficiency and user satisfaction



CASSETTE TYPE
Features a compact design, allowing users to collect a urine sample and apply it to the testing area

PREGNANCY TEST KITS

MIDSTREAM TYPE



Diagnostic
Solutions



FEATURES



HOME
USE



MIDSTREAM
TYPE



QUICK &
ACCURATE



BEST
SELLER

AME Smart Stream Pregnancy Test Kit Mid Stream Type I

DETAILED FEATURES

- Results in 40 seconds Highest in class test area
- Quick and accurate results with accuracy of >99.9%
- Bagged in Aluminium Jacket
- Includes Dual Language Leaflet
- Home Use
- Midstream Type



Diagnostic Solutions

PREGNANCY TEST KITS CASSETTE & MIDSTREAM TYPE



AME SMART CHECK

AME DIGITAL SMART STREAM

TECHNOLOGY



FEATURES

Quick and accurate results with accuracy of >99.9%
Suitable for Home Use

Results in 5 Minutes Highest in class test area
Cassette type
Bagged in Aluminium Jacket
Includes Dual Language Leaflet

Self-testing format for High Privacy
Convenient "Yes" / "No" Indicator
Easy to use with 1 Step Operation
Over 99.9% Accuracy
Fast Results in 3 Minutes
Simple visual interpretation
Shelf Life 2 Years
Sensitivity >99.9%
Specificity >99.9%
Cut-Off 25 mIU/mL

FERTILITY TEST KITS

MIDSTREAM TYPE



Diagnostic Solutions



AME SMART OVU TEST



AME SMART STREAM MENOPAUSE

TECHNOLOGY



Home USE



Quick & ACCURATE



Midstream STYLE

FEATURES

Quick and accurate results with accuracy of >99.9%
Suitable for Home Use

AME SmartOvuTest enables couples trying to conceive for a baby to try Smartly. Utilizing the Menstrual Cycle, the test helps couples find the most fertile days.

5 Test kits in 1 pack.

Suitable for daily testing

Easy to use

For any cycle length. One pack is sufficient for most women with regular cycles.

Additional tests maybe needed for irregular cycles

Discover your likely menopausal stage - By combining FSH test results, cycle history, and symptoms, the test can indicate if a woman is likely to have Menopause or Not. Key Features Include:

Easy to Interpret Results.

Simple, Reliable, and Accurate with FSH Hormone sensitivity of 25 mIU/mL.

Rapid & Simple to Use - Results in 5 Minutes

Easy to Read Results - two lines highlight a positive result, just like your regular HCG test kits.

DRUG OF ABUSE TESTS

Welcome to the forefront of personalized diagnostics with AME's In-Vitro Diagnostics Kits. Our Drug of Abuse Test Kits redefine convenience and accuracy in home testing. Designed with precision, AME empowers users with reliable results, enabling proactive health management. Experience a new era in at-home healthcare with AME's In-Vitro Diagnostics Kits - where accuracy meets convenience.



FEATURES

Our products are designed with a focus on accuracy, safety, and user convenience, ensuring reliable performance. Each device is manufactured to meet international quality and regulatory standards, providing healthcare professionals and consumers with confidence in every use.



HOME USE

Provides convenience and accurate results in the comfort of your home



CASSETTE TYPE

Features a compact design, allowing users to collect a urine sample and apply it to the testing area



QUICK & ACCURATE

Ensures reliable results in minimal time, enhancing efficiency and user satisfaction

DRUG OF ABUSE TESTS

MIDSTREAM TYPE



Diagnostic Solutions



AME MULTI DRUG OF ABUSE KIT

4 in 1 Test Kit (Cassette Type)

Saliva Test sample

Designed to test four common drugs of abuse



Cassette
STYLE



Quick &
ACCURATE



AME MULTI DRUG OF ABUSE KIT

6 in 1 Test Kit (Cassette Type)

Urine Test sample

Designed to test six common drugs of abuse



Cassette
STYLE



Quick &
ACCURATE



AME MULTI DRUG OF ABUSE KIT

Cup Type

6 to 15 Parameters



Cassette
STYLE



Quick &
ACCURATE



EASY
TO USE



QUICK AND
ACCURATE



OVER 99%
ACCURACY



AME ALCOHOL RAPID TEST

Saliva sample

Easy to use

Quick and accurate



Cassette
STYLE



Quick &
ACCURATE

Maternity Care Solutions



Our Maternity Care Solutions are designed to support women through every stage of their reproductive journey, from planning conception to confirming pregnancy with confidence. Combining precision diagnostics with user-friendly design, our products deliver fast, reliable results for both clinical and home use. We aim to empower women and healthcare providers alike with dependable tools that promote informed decisions and better care outcomes.

Product Range



Breast Pumps



BREAST PUMPS

Introducing the Breast Pumps category, designed to enhance the breastfeeding experience for modern mothers. Our state-of-the-art pumps combine comfort, efficiency, and portability, providing a seamless solution for nursing mothers on the go. With customizable settings and advanced technology, trust in AME's commitment to nurturing maternal well-being, as our innovative design ensures a discreet and efficient pumping experience.

CATEGORIES

-  Electric Breast Pump
-  Manual Breast Pump



FEATURES

Our products are designed with a focus on accuracy, safety, and user convenience, ensuring reliable performance. Each device is manufactured to meet international quality and regulatory standards, providing healthcare professionals and consumers with confidence in every use.



Small sized device with multiple functions; 4 modes and 9 levels



Silent & discreet pumping



Stops the milks back flow in the device



Compact and lightweight



Safe to use material



A soft LED light, for night time pumping



Easy to Assemble



Different colors available



easy to clean the device and parts



Easy to clean spills



Good suction function



Large suction capacity



Saves time with quick and efficient pumping



Maternity Solutions

BREAST PUMPS ELECTRIC & MANUAL



AME LD 208
Electric Model
Wearable Type

FEATURES



Small Size
Big Functions,
4 Modes, 9 Levels



Easy to
Assemble



Take it
Everywhere



Colorful
Design

DETAILED FEATURES

- Wearable design
- 6 levels of stimulate modes
- 9 levels of suction modes
- 9 levels of expression modes
- 9 levels of multi suction modes
- Breast shield (flange size) 24/27mm
- Silent and discreet pumping
- Compact and lightweight
- Powered by Lithium 1200mAh/TYPE-C 5V, 1A (include cable)

FEATURES



Easy to
Clean Spills



Easy to
Clean



BPA
Free



Large
Suction



Save
Time

DETAILED FEATURES

- Streamline and comfortable handle
- Soft and comfortable breast-shield to fit snugly against breast
- Silent and discreet pumping
- Easy to assemble, use and clean
- Unique system for easy milk storage
- Compact and lightweight
- BPA free



AME LD 101
Manual Model

BREAST PUMPS ELECTRIC



Maternity
Solutions



AME LD 202

AME LD 209

TECHNOLOGY



Easy to
Clean



Easy to
Assemble



BPA
Free



Small Size
Big Functions



Ultra-
Quiet



Anti-
BACKFLOW



Colorful
Design



Reasonable
suction

FEATURES

Silent and discreet pumping
Compact and lightweight
BPA free

2-phase expression technology
Powerful motor, 10 suction level options
Soft and comfortable breast-shield
Easy to assemble, use and clean
Unique system for easy milk storage
Powered by 4AA batteries or adaptor 6V/1A

6 stimulate modes; 9 suction modes;
9 expression modes; 9 multi suction modes
24 mm comfortable breastshield
Automatically power-off in 30 minutes
Powered by 5V/1A type-c adaptor/cable
(only include cable)



Maternity
Solutions

BREAST PUMPS ELECTRIC



AME LD 302

AME LD 305

TECHNOLOGY



Easy to
Clean



Easy to
Assemble



BPA
Free



Small Size
Big Functions



Ultra-
Quiet

FEATURES

- 2-phase expression technology
- Powerful motor, 10 suction level options
- Dual Cup type
- Soft and comfortable breast-shield to fit snugly against breast
- Silent and discreet pumping
- Easy to assemble, use and clean
- Unique system for easy milk storage
- Compact and lightweight
- BPA free
- Powered by 4AA batteries or adaptor 6V/1A

BREAST PUMPS

ELECTRIC



Maternity
Solutions



AME LD 3010L

AME LD 2010

TECHNOLOGY



Colorful
Design



Reasonable
suction



LED
Night Lamp

FEATURES

Dual Cup type
 Silent and discreet pumping
 Compact and lightweight
 BPA free
 6 stimulate modes; 9 suction modes; 9 expression modes; 9 multi suction modes
 24 mm comfortable breast-shield
 Automatically power-off in 30 minutes

Powered by Li-ion 2000mAh/TYPE-C 5V, 1A
(Rechargeable Type)
 LED lamp (Night Light)

Powered by 5V/1A type-c adaptor/cable/
 4"AA"batteries (only include cable)

Pain Relief & Wellness Solutions



Our Pain Relief & Wellness Solutions are crafted to enhance comfort, recovery, and overall well-being. From soothing pain patches to revitalizing wellness aids, each product is designed with innovative materials and technology to deliver effective, non-invasive relief. Whether for clinical use or everyday comfort, our solutions support a healthier, more balanced lifestyle.

Product Range



PAIN & COOLING PATCHES

Discover the versatility and reliability of AME's Medical and Non-Medical Patches, designed to cater to a wide range of healthcare and personal needs. Our medical patches provide targeted therapeutic solutions, ensuring precise drug delivery and enhanced patient care. For non-medical applications, AME offers innovative patches that prioritize comfort and durability for everyday use. Engineered with high-quality materials, our patches are skin-friendly, adhesive, and built for performance. Trust AME for cutting-edge patch solutions that combine science and innovation to meet diverse requirements effectively.



FEATURES

Our products are designed with a focus on accuracy, safety, and user convenience, ensuring reliable performance. Each device is manufactured to meet international quality and regulatory standards, providing healthcare professionals and consumers with confidence in every use.



COOLING PATCH

Provide instant, drug-free relief by reducing skin temperature and offering soothing comfort



PAIN RELIEF PATCH

Deliver targeted, long-lasting relief by directly soothing sore muscles and joints



CONTINUOUS HEATING

Provide continuous warmth to relax muscles and alleviate discomfort for hours



COLOR TRANSITION

Color transition provides a visual indicator of fever from blue at normal temperature to pink at higher body temperatures



Wellness
Solutions

PAIN & COOLING PATCHES

COOLING PATCH



FEATURES



COOLING
PATCH



COLOR
TRANSITION

AME MEDICAL COOLING PATCH Child

DETAILED FEATURES

AME Medical Cooling gel patch with antipyretic cooling effect; can relieve fever, headache and toothache.

Thermal discoloration - Blue to Red

Super viscosity

Non-allergic

Good cooling effect and Sticking Tendency.

Menthol Level customization as per customer requirements.

Can be produced for Adult & Children alike.

Pack Sizes Available: 4's, 6's, 25's and 50's



AME MEDICAL COOLING PATCH Adult

PAIN & COOLING PATCHES

PAIN PATCH



Wellness Solutions



PAIN RELIEF Patch



CONTINUOUS Heating



AME FAR INFRARED PATCH

Experience the ancient wisdom of herbal remedies combined with cutting-edge Far Infrared Technology.

Our Pain Relief Heat Patch harnesses the power of nature to soothe your aches and pains effortlessly. Infused with a blend of carefully selected herbal extracts, our patches offer a holistic approach to pain management. Let the natural goodness of herbs work its magic on your body, promoting healing from within.

Sleek 8cm x 10cm size – perfect for on-the-go relief. Pack Sizes Available: 4's, 6's, 25's and 50's



PAIN RELIEF Patch



CONTINUOUS Heating



AME MENSTRUAL CRAMP PATCHES

AME Menstrual Cramp Pain relief Patches enable deep heat into the pelvic muscles and accelerate blood circulation to provide fast and effective comfort for the lower abdomen.

Product Benefits:

Continuous heating, for up to 8 hours, providing long-lasting warmth and warmth.

Made from Natural Herbal Extracts - Drug Free



PAIN RELIEF Patch



CONTINUOUS Heating



AME CAPSICUM HEAT PATCHES

AME Capsicum Heat Patches can relieve pain symptoms caused by bruising, arthritis, Rheumatism, general sprains and strain and muscle related back pain.

Product Benefits:

Capsaicin, the primary compound in hot peppers and proven to reduce discomfort.

12 Hour Continuous Use

Available in pack of 6 patches

8 x 13cm

Suitable for adults and children over 12 years

FOOT CREAM

Introducing AME's Foot Cream, a therapeutic blend enriched with Urea, Vitamin E, AHA's, and Dimethicone, meticulously crafted to pamper and rejuvenate tired feet. This unique formulation offers deep moisturization, promoting softness and suppleness while addressing dryness and roughness. AME's Foot Cream stands as a testament to our commitment to foot care excellence, providing a nourishing solution that combines science-backed ingredients for optimal skin health. Trust in AME for a soothing experience that goes beyond the surface, ensuring your feet receive the care they deserve. Step into comfort with AME Foot Cream – where indulgence meets efficacy.



AME FOOT CREAM 10% UREA

FEATURES

AME's Foot Cream is a therapeutic blend enriched with:
10% Urea for Moisture Retention
Dimethicone as Moisturizing Agent
AHA's to remove dead skin
Vitamin E for Glowing Skin

It is meticulously crafted to pamper and rejuvenate tired feet. This unique formulation offers deep moisturization, promoting softness and suppleness while addressing dryness and roughness.

AME's Foot Cream stands as a testament to our commitment to foot care excellence, providing a nourishing solution that combines science-backed ingredients for optimal skin health.



Safe Injections Solutions



Ensure safety and precision in healthcare with AME's Auto Disable Syringes, designed to prevent reuse and minimize infection risks. Engineered with cutting-edge technology, our syringes automatically lock after a single use, ensuring maximum patient and healthcare worker protection. AME's Auto Disable Syringes combine reliability with ease of use, meeting the highest safety and quality standards. Ideal for vaccination programs and routine medical procedures, they are a trusted choice for safe and efficient care. Choose AME for innovative solutions that prioritize safety without compromise.



STEP I DRAW

If using the unmounted version, make sure the needle is fixed tight and then use the standard aseptic technique to DRAW out the medication and fill the syringe.



STEP II INJECT

INJECT the medication by fully depressing the plunger to the end of the syringe barrel to engage the locking mechanism. The user should be able to feel and hear a 'CLICK'.



STEP III RETRACT

ZIP the plunger back to retract the needle safely back into the barrel. The needle will tilt to one side and encapsulate the needle in the barrel of the syringe, preventing reuse or needlestick injury.



STEP IV BREAK

SNAP off the plunger at the breaking point. Your AME AD syringe is now disabled and cannot be reused. Discard the complete unit as per regulations require

AUTO DISABLE SYRINGE SMARTSTEPS™ SYSTEM



Safe Injection
Solutions

TECHNOLOGY



**MODEL: AME STANDARD
AUTO DISABLE SYRINGE**

FEATURES

23-31G, needle length 6mm-25mm, thin wall/regular wall
Sterile, non-toxic, non-pyrogenic
Single use only
Shelf life: 5 years

CERTIFICATION

Certification and Quality Assurance: MDR and FDA 510k approved and manufactured in accordance with ISO 13485

SPECIFICATIONS

Main Material: PE, PP, PC, SUS304 Stainless Steel
Cannula, Silicone Oil
Material for gasket: Isoprene rubber, Latex free
Safety design and easy for using



**MODEL: AME INSULIN
SYRINGE**

FEATURES

Structure and Composition: The Syringes are assembled by sliding sleeve, needle protective cap, needle tube, barrel, plunger and piston. All parts and material for this product meet medical requirements. Sterilized by EO. Intend Use: The sterile syringes for insulin for single use intended to be used to administer an injection of insulin to a patient subcutaneously.

CERTIFICATION

Certificate and Quality Assurance: CE, FDA approved and manufactured in accordance with ISO 13485

SPECIFICATIONS

Main Material: PP, isoprene rubber, silicone oil, SUS304 stainless steel cannula.
Specification: 1ml, 0.5ml, 0.3ml, U-40, U-100.
Needle Size: 27G-31G.

Surgical Tools & Hospital Furniture Solutions



Our Surgical Instruments & Hospital Furniture Solutions are built to deliver precision, durability, and comfort in every healthcare setting. From high-quality surgical tools to ergonomically designed hospital furniture, our range supports efficient workflows and patient-centered care. Each product is crafted to meet rigorous medical standards, ensuring reliability and performance in every procedure.

Product Range



**Surgical
Instruments**



**Hospital
Furniture**



HOSPITAL FURNITURE

Step into a world of comfort and functionality with AME's Hospital Furniture, meticulously designed to enhance patient care and optimize healthcare environments. From durable patient beds to ergonomic seating solutions, our furniture seamlessly blends reliability with contemporary design. AME's commitment to excellence ensures that each piece is crafted for durability, infection control, & patient well-being. Trust in AME to provide hospital furniture that aligns with the evolving needs of modern healthcare settings, creating spaces that prioritize both comfort and efficiency. Elevate your healthcare infrastructure with AME.

PRODUCT LINE

- Hospital Beds
- Wheelchairs
- Bathroom Safety Products
- Mobility Aid Products



Please download this range by scanning the QR Code

SURGICAL INSTRUMENTS

Embark on precision and excellence in surgery with AME's Surgical Instruments. Crafted with the highest quality materials and precision engineering, our instruments ensure optimal performance and reliability in the operating room. From scalpels to forceps, AME's comprehensive range caters to the diverse needs of surgical procedures. Trust in the expertise of AME to deliver instruments that embody the highest standards of safety, precision, and durability, fostering confidence among healthcare professionals. Join us in advancing surgical precision with AME – where every instrument reflects a commitment to exceptional surgical outcomes.

PRODUCT LINE

- Operating Scissors
- Needle Holders
- Dressing Forceps
- Stainless steel Hollowware
- Hemostatic Forceps
- Diagnostic Instruments



Please download this range by scanning the QR Code

CONTACT US

**AT AME GROUP
HEAD OFFICE**



Reach Us

+971 6 558 4052

connect@ame.ae

Visit Our Website

www.ameworldwide.com

Visit Our Office

Warehouse no. 6, Street 47,
Industrial Area 13, Sharjah, UAE

www.almaniads.com

